

Antigen peptide transporter 2

Catalog No: #AP88055



Package Size: #AP88055-1 50ug #AP88055-2 100ug #AP88055-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Antigen peptide transporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	APT2,ATP-binding cassette sub-family B member 3,Histocompatibility antigen modifier 2,Abcb3,Ham-2,Ham2
Accession No.	Uniprot ID:P36371
Uniprot	P36371
GeneID	21355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tap2

Note: This product is for in vitro research use only