

# Selenocysteine insertion sequence-binding protein 2-like

Catalog No: #AP88059

Package Size: #AP88059-1 50ug #AP88059-2 100ug #AP88059-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Selenocysteine insertion sequence-binding protein 2-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SECIS-binding protein 2-like,Kiaa0256
Accession No.	Uniprot ID:Q6A098
Uniprot	Q6A098
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Secisbp2l

Note: This product is for in vitro research use only