

HEN2 Antibody

Catalog No: #44478

Package Size: #44478-1 50ul #44478-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HEN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	HEN2 Antibody detects endogenous levels of total HEN2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HEN2
Target Name	HEN2
Other Names	NHLH2, BHLHa34, HEN-2, KIAA0490, Helix-loop-helix protein 2, HEN2, Nescient helix loop helix 2, NSCL-2, NSCL2
Accession No.	Swiss-Prot#: Q02577NCBI Gene ID:4808
Uniprot	Q02577
GeneID	4808;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only