

Sperm surface protein Sp17

Catalog No: #AP88060



Package Size: #AP88060-1 50ug #AP88060-2 100ug #AP88060-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sperm surface protein Sp17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Sperm autoantigenic protein 17,Sp17
Accession No.	Uniprot ID:Q62252
Uniprot	Q62252
GeneID	20686;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spa17

Note: This product is for in vitro research use only