

Reticulon-4

Catalog No: #AP88065



Package Size: #AP88065-1 50ug #AP88065-2 100ug #AP88065-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Reticulon-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurite outgrowth inhibitor, Nogo protein, Kiaa0886, Nogo
Accession No.	Uniprot ID:Q99P72
Uniprot	Q99P72
GeneID	68585;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rtn4

Note: This product is for in vitro research use only