

Sphingosine 1-phosphate receptor 1

Catalog No: #AP88068



Package Size: #AP88068-1 50ug #AP88068-2 100ug #AP88068-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	S1P receptor 1,Endothelial differentiation G-protein coupled receptor 1,Sphingosine 1-phosphate receptor Edg-1,S1P receptor Edg-1,CD363,Edg1
Accession No.	Uniprot ID:P48303
Uniprot	P48303
GeneID	29733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S1pr1

Note: This product is for in vitro research use only