

Retinoic acid receptor responder protein 3

Catalog No: #AP88069



Package Size: #AP88069-1 50ug #AP88069-2 100ug #AP88069-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor responder protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HRAS-like suppressor 4,HRSL4,RAR-responsive protein TIG3,Retinoid-inducible gene 1 protein,Tazarotene-induced gene 3 protein,RIG1,TIG3
Accession No.	Uniprot ID:Q9UL19
Uniprot	Q9UL19
GeneID	5920;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RARRES3

Note: This product is for in vitro research use only