

Somatostatin receptor type 5

Catalog No: #AP88074



Package Size: #AP88074-1 50ug #AP88074-2 100ug #AP88074-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SS-5-R,Smstr5
Accession No.	Uniprot ID:O08858
Uniprot	O08858
GeneID	20609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sstr5

Note: This product is for in vitro research use only