

Serine/arginine-rich splicing factor 3

Catalog No: #AP88076



Package Size: #AP88076-1 50ug #AP88076-2 100ug #AP88076-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Splicing factor, arginine/serine-rich 3,SFRS3
Accession No.	Uniprot ID:Q3SZR8
Uniprot	Q3SZR8
GeneID	540276;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRSF3

Note: This product is for in vitro research use only