

# Serine/threonine-protein kinase 17B

Catalog No: #AP88077



Package Size: #AP88077-1 50ug #AP88077-2 100ug #AP88077-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase 17B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DAP kinase-related apoptosis-inducing protein kinase 2,Drak2
Accession No.	Uniprot ID:Q91XS8
Uniprot	Q91XS8
GeneID	170904;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Stk17b

Note: This product is for in vitro research use only