

Sulfotransferase 1A1

Catalog No: #AP88078



Package Size: #AP88078-1 50ug #AP88078-2 100ug #AP88078-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sulfotransferase 1A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ST1A1,Aryl sulfotransferase,Phenol sulfotransferase,Phenol/aryl sulfotransferase,mSTp1,ST1A4,Sulfokinase,St1a1,Stp,Stp1
Accession No.	Uniprot ID:P52840
Uniprot	P52840
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sult1a1

Note: This product is for in vitro research use only