

Semaphorin-4D

Catalog No: #AP88081



Package Size: #AP88081-1 50ug #AP88081-2 100ug #AP88081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin-4D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	A8,BB18,GR3,CD100,C9orf164,CD100,SEMAJ
Accession No.	Uniprot ID:Q92854
Uniprot	Q92854
GeneID	10507;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SEMA4D

Note: This product is for in vitro research use only