

TATA-box-binding protein

Catalog No: #AP88086



Package Size: #AP88086-1 50ug #AP88086-2 100ug #AP88086-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TATA-box-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TATA sequence-binding protein,TATA-binding factor,TATA-box factor,Transcription initiation factor TFIID TBP subunit
Accession No.	Uniprot ID:Q2HJ52
Uniprot	Q2HJ52
GeneID	516578;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TBP

Note: This product is for in vitro research use only