

# TGF-beta-activated kinase 1 and MAP3K7-binding protein 3

Catalog No: #AP88089

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP88089-1 50ug #AP88089-2 100ug #AP88089-3 1mg

## Description

|                    |   |
|--------------------|---|
| Product Name       | TGF-beta-activated kinase 1 and MAP3K7-binding protein 3  |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | Mitogen-activated protein kinase kinase kinase 7-interacting protein 3,NF-kappa-B-activating protein 1,TAK1-binding protein 3,TAB-3,TGF-beta-activated kinase 1-binding protein 3,MAP3K7IP3 |
| Accession No.      | Uniprot ID:Q8N5C8   |
| Uniprot            | Q8N5C8  |
| GeneID             | 257397;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TAB3

Note: This product is for in vitro research use only