

Suppressor of cytokine signaling 2

Catalog No: #AP88091



Package Size: #AP88091-1 50ug #AP88091-2 100ug #AP88091-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	SOCS-2
Accession No.	Uniprot ID:Q7YRV6
Uniprot	Q7YRV6
GeneID	100037966;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SOCS2

Note: This product is for in vitro research use only