

Secreted frizzled-related protein 1

Catalog No: #AP88095



Package Size: #AP88095-1 50ug #AP88095-2 100ug #AP88095-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secreted frizzled-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FRP-1,Secreted apoptosis-related protein 2,SARP-2,FRP,FRP1,SARP2
Accession No.	Uniprot ID:Q8N474
Uniprot	Q8N474
GeneID	6422;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFRP1

Note: This product is for in vitro research use only