

Sodium-dependent serotonin transporter

Catalog No: #AP88105



Package Size: #AP88105-1 50ug #AP88105-2 100ug #AP88105-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium-dependent serotonin transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SERT,5HT transporter,5HTT,Solute carrier family 6 member 4,Htt,Sert
Accession No.	Uniprot ID:Q60857
Uniprot	Q60857
GeneID	15567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc6a4

Note: This product is for in vitro research use only