

Switch-associated protein 70

Catalog No: #AP88106

Package Size: #AP88106-1 50ug #AP88106-2 100ug #AP88106-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Switch-associated protein 70
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SWAP-70,Kiaa0640
Accession No.	Uniprot ID:Q6A028
Uniprot	Q6A028
GeneID	20947;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Swap70

Note: This product is for in vitro research use only