

Lysine--tRNA ligase

Catalog No: #AP88107



Package Size: #AP88107-1 50ug #AP88107-2 100ug #AP88107-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysine--tRNA ligase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lysyl-tRNA synthetase,LysRS,KIAA0070
Accession No.	Uniprot ID:Q15046
Uniprot	Q15046
GeneID	3735;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KARS

Note: This product is for in vitro research use only