

Somatostatin receptor type 2

Catalog No: #AP88111



Package Size: #AP88111-1 50ug #AP88111-2 100ug #AP88111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SS-2-R,SRIF-1,Smstr2,Sst2
Accession No.	Uniprot ID:P30875
Uniprot	P30875
GeneID	20606;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sstr2

Note: This product is for in vitro research use only