

T-box transcription factor TBX1

Catalog No: #AP88113



Package Size: #AP88113-1 50ug #AP88113-2 100ug #AP88113-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-box transcription factor TBX1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	T-box protein 1, Testis-specific T-box protein
Accession No.	Uniprot ID:O43435
Uniprot	O43435
GeneID	6899;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TBX1

Note: This product is for in vitro research use only