

Protein Spindly

Catalog No: #AP88114



Package Size: #AP88114-1 50ug #AP88114-2 100ug #AP88114-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein Spindly
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Coiled-coil domain-containing protein 99,Spindle apparatus coiled-coil domain-containing protein 1,Ccdc99
Accession No.	Uniprot ID:Q3KR99
Uniprot	Q3KR99
GeneID	303037;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spd11

Note: This product is for in vitro research use only