

Signal transducer and activator of transcription 5A

Catalog No: #AP88115



Package Size: #AP88115-1 50ug #AP88115-2 100ug #AP88115-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Mammary gland factor, MGF
Accession No.	Uniprot ID:Q95115
Uniprot	Q95115
GeneID	282375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: STAT5A

Note: This product is for in vitro research use only