

Reticulon-4

Catalog No: #AP88119



Package Size: #AP88119-1 50ug #AP88119-2 100ug #AP88119-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Reticulon-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Foocen, Glut4 vesicle 20 kDa protein, Neurite outgrowth inhibitor, Nogo protein, Nogo
Accession No.	Uniprot ID:Q9JK11
Uniprot	Q9JK11
GeneID	83765;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rtn4

Note: This product is for in vitro research use only