

Metalloproteinase inhibitor 1

Catalog No: #AP88120



Package Size: #AP88120-1 50ug #AP88120-2 100ug #AP88120-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Tissue inhibitor of metalloproteinases 1, TIMP-1
Accession No.	Uniprot ID:P81546
Uniprot	P81546
GeneID	403816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TIMP1

Note: This product is for in vitro research use only