

Tyrosine-protein phosphatase non-receptor type substrate 1

Catalog No: #AP88124

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88124-1 50ug #AP88124-2 100ug #AP88124-3 1mg

Description

Product Name	Tyrosine-protein phosphatase non-receptor type substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SHP substrate 1,Brain Ig-like molecule with tyrosine-based activation motifs,Bit,CD172 antigen-like family member A,Inhibitory receptor SHPS-1,Macrophage fusion receptor,Macrophage membrane protein MFP150,Signal-regulatory protein alpha-1,Sirp-alpha-1,CD172a,Bit,Mfr,Ptpns1,Shps1,Sirp
Accession No.	Uniprot ID:P97710
Uniprot	P97710
GeneID	25528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sirpa

Note: This product is for in vitro research use only