

Antileukoproteinase

Catalog No: #AP88126



Package Size: #AP88126-1 50ug #AP88126-2 100ug #AP88126-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Antileukoproteinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ALP,BLPI,HUSI-1,Mucus proteinase inhibitor,MPI,Protease inhibitor WAP4,Secretory leukocyte protease inhibitor,Seminal proteinase inhibitor,WAP four-disulfide core domain protein 4,WAP4,WFDC4
Accession No.	Uniprot ID:P03973
Uniprot	P03973
GeneID	6590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLPI

Note: This product is for in vitro research use only