

Toll-like receptor 4

Catalog No: #AP88127



Package Size: #AP88127-1 50ug #AP88127-2 100ug #AP88127-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toll-like receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Toll4,CD284
Accession No.	Uniprot ID:Q9QX05
Uniprot	Q9QX05
GeneID	29260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tlr4

Note: This product is for in vitro research use only