

O-phosphoseryl-tRNA(Sec) selenium transferase

Catalog No: #AP88129



Package Size: #AP88129-1 50ug #AP88129-2 100ug #AP88129-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	O-phosphoseryl-tRNA(Sec) selenium transferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Selenocysteine synthase, Sec synthase, SelenocysteinyI-tRNA(Sec) synthase, Sep-tRNA:Sec-tRNA synthase, SepSecS, UGA suppressor tRNA-associated protein, D5ErtD135e
Accession No.	Uniprot ID:Q6P6M7
Uniprot	Q6P6M7
GeneID	211006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sepsecs

Note: This product is for in vitro research use only