

Sestrin-1

Catalog No: #AP88130



Package Size: #AP88130-1 50ug #AP88130-2 100ug #AP88130-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sestrin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pa26,Sest1
Accession No.	Uniprot ID:P58006
Uniprot	P58006
GeneID	140742;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sesn1

Note: This product is for in vitro research use only