

## Serine/arginine repetitive matrix protein 4

Catalog No: #AP88131



Package Size: #AP88131-1 50ug #AP88131-2 100ug #AP88131-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/arginine repetitive matrix protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Medulloblastoma antigen MU-MB-2.76,Neural-specific serine/arginine repetitive splicing factor of 100 kDa,Neural-specific SR-related protein of 100 kDa,KIAA1853
Accession No.	Uniprot ID:A7MD48
Uniprot	A7MD48
GeneID	84530;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SRRM4

Note: This product is for in vitro research use only