

Superoxide dismutase [Cu-Zn]

Catalog No: #AP88133



Package Size: #AP88133-1 50ug #AP88133-2 100ug #AP88133-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Superoxide dismutase [Cu-Zn]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:P04178
Uniprot	P04178
GeneID	397036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SOD1

Note: This product is for in vitro research use only