

## Stabilin-2

Catalog No: #AP88138



Package Size: #AP88138-1 50ug #AP88138-2 100ug #AP88138-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Stabilin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fasciclin, EGF-like, laminin-type EGF-like and link domain-containing scavenger receptor 2,FEEL-2,Feel2,Hare
Accession No.	Uniprot ID:Q8R4U0
Uniprot	Q8R4U0
GeneID	192188;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Stab2

Note: This product is for in vitro research use only