

TGF-beta receptor type-1

Catalog No: #AP88142



Package Size: #AP88142-1 50ug #AP88142-2 100ug #AP88142-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF-beta receptor type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TGFR-1,TGF-beta type I receptor,Transforming growth factor-beta receptor type I,TGF-beta receptor type I
Accession No.	Uniprot ID:O46680
Uniprot	O46680
GeneID	282382;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFBR1

Note: This product is for in vitro research use only