

THAP domain-containing protein 4

Catalog No: #AP88143



Package Size: #AP88143-1 50ug #AP88143-2 100ug #AP88143-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	THAP domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q6P3Z3
Uniprot	Q6P3Z3
GeneID	67026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Thap4

Note: This product is for in vitro research use only