

Protein S100-A9

Catalog No: #AP88146



Package Size: #AP88146-1 50ug #AP88146-2 100ug #AP88146-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Protein S100-A9 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Calgranulin-B, Migration inhibitory factor-related protein 14, MRP-14, Myeloid-related protein 14, S100 calcium-binding protein A9, Mrp14 |
| Accession No. | Uniprot ID:P50116 |
| Uniprot | P50116 |
| GeneID | 94195; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100a9

Note: This product is for in vitro research use only