

Tenascin-R

Catalog No: #AP88149



Package Size: #AP88149-1 50ug #AP88149-2 100ug #AP88149-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenascin-R
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TN-R,Janusin,Neural recognition molecule J1-160/180,Restrictin
Accession No.	Uniprot ID:Q8BYI9
Uniprot	Q8BYI9
GeneID	21960;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnr

Note: This product is for in vitro research use only