

Sphingosine 1-phosphate receptor 3

Catalog No: #AP88152

Package Size: #AP88152-1 50ug #AP88152-2 100ug #AP88152-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	S1P receptor 3,Endothelial differentiation G-protein coupled receptor 3,Lysophospholipid receptor B3,Sphingosine 1-phosphate receptor Edg-3,S1P receptor Edg-3,Edg3,Lpb3
Accession No.	Uniprot ID:Q9Z0U9
Uniprot	Q9Z0U9
GeneID	13610;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S1pr3

Note: This product is for in vitro research use only