

Secretoglobin family 2A member 2

Catalog No: #AP88154



Package Size: #AP88154-1 50ug #AP88154-2 100ug #AP88154-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secretoglobin family 2A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C3.1,Prostatein peptide C3,Prostatic steroid-binding protein C3,Psbpc3
Accession No.	Uniprot ID:P02780
Uniprot	P02780
GeneID	25010;361725;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Scgb2a2

Note: This product is for in vitro research use only