

Ras-related protein Rab-8A

Catalog No: #AP88155



Package Size: #AP88155-1 50ug #AP88155-2 100ug #AP88155-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rab-8A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_2k8
Accession No.	Uniprot ID:Q5F470
Uniprot	Q5F470
GeneID	428352;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RAB8A

Note: This product is for in vitro research use only