

GTP-binding nuclear protein Ran

Catalog No: #AP88156



Package Size: #AP88156-1 50ug #AP88156-2 100ug #AP88156-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GTP-binding nuclear protein Ran
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	GTPase Ran,Ras-like protein TC4,Ras-related nuclear protein
Accession No.	Uniprot ID:P62825
Uniprot	P62825
GeneID	100499482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAN

Note: This product is for in vitro research use only