

T-cell leukemia/lymphoma protein 1A

Catalog No: #AP88158



Package Size: #AP88158-1 50ug #AP88158-2 100ug #AP88158-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell leukemia/lymphoma protein 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Oncogene TCL-1,Oncogene TCL1,Protein p14 TCL1,TCL1
Accession No.	Uniprot ID:P56279
Uniprot	P56279
GeneID	8115;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TCL1A

Note: This product is for in vitro research use only