

Secreted frizzled-related protein 1

Catalog No: #AP88163



Package Size: #AP88163-1 50ug #AP88163-2 100ug #AP88163-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secreted frizzled-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	sFRP-1,Frizzled in aorta protein,FrzA protein,FRZA
Accession No.	Uniprot ID:O19116
Uniprot	O19116
GeneID	282068;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFRP1

Note: This product is for in vitro research use only