

Trefoil factor 2

Catalog No: #AP88166



Package Size: #AP88166-1 50ug #AP88166-2 100ug #AP88166-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trefoil factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Pancreatic spasmolytic polypeptide,PSP,Spasmolytic polypeptide,SP,SML1
Accession No.	Uniprot ID:P01359
Uniprot	P01359
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TFF2

Note: This product is for in vitro research use only