

Transcription initiation factor TFIID subunit 1

Catalog No: #AP88167



Package Size: #AP88167-1 50ug #AP88167-2 100ug #AP88167-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription initiation factor TFIID subunit 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell cycle gene 1 protein, TBP-associated factor 250 kDa, p250, Transcription initiation factor TFIID 250 kDa subunit, TAF(II)250, BA2R, CCG1, CCGS, TAF2A
Accession No.	Uniprot ID: P21675
Uniprot	P21675
GeneID	6872;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TAF1

Note: This product is for in vitro research use only