

Serpín B5

Catalog No: #AP88169



Package Size: #AP88169-1 50ug #AP88169-2 100ug #AP88169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serpín B5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Maspin,Peptidase inhibitor 5,PI-5,PI5
Accession No.	Uniprot ID:P70564
Uniprot	P70564
GeneID	116589;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpínb5

Note: This product is for in vitro research use only