

Transforming growth factor beta-2

Catalog No: #AP88176



Package Size: #AP88176-1 50ug #AP88176-2 100ug #AP88176-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TGF-beta-2
Accession No.	Uniprot ID:Q07257
Uniprot	Q07257
GeneID	81809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tgfb2

Note: This product is for in vitro research use only