

Thyroid hormone receptor beta

Catalog No: #AP88178



Package Size: #AP88178-1 50ug #AP88178-2 100ug #AP88178-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyroid hormone receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Nuclear receptor subfamily 1 group A member 2,C-ERBA-BETA,NR1A2
Accession No.	Uniprot ID:P68306
Uniprot	P68306
GeneID	396431;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:THRB

Note: This product is for in vitro research use only