

Metalloproteinase inhibitor 2

Catalog No: #AP88179



Package Size: #AP88179-1 50ug #AP88179-2 100ug #AP88179-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Tissue inhibitor of metalloproteinases 2,TIMP-2
Accession No.	Uniprot ID:Q9TTY1
Uniprot	Q9TTY1
GeneID	403633;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TIMP2

Note: This product is for in vitro research use only