

Transaldolase

Catalog No: #AP88181



Package Size: #AP88181-1 50ug #AP88181-2 100ug #AP88181-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transaldolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TAL,TALDO,TALDOR
Accession No.	Uniprot ID:P37837
Uniprot	P37837
GeneID	6888;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TALDO1

Note: This product is for in vitro research use only